

## SECULIFE | VL Variable Load

3-349-872-03  
1/7.15

### Functional testing of Defibrillators

- Change the variable load of SECULIFE | DF<sub>PRO/BASE</sub>
- Capture of energy output
- Fully AED Compatible
- Easy to use



### Applications

The SECULIFE | VL makes it possible to adjust the variable loads of the SECULIFE | DF<sub>PRO/BASE</sub> from 25 Ω to 200 Ω and back down again in 25 Ω steps by simply turning a knob.

Connection via the RS 232 port permits automatic detection of the selected external load.

### Overall structure



### Variable-Load-Table

25 Ω → 50 Ω
50 Ω → 75 Ω
75 Ω → 100 Ω
100 Ω → 125 Ω
125 Ω → 150 Ω
150 Ω → 175 Ω
175 Ω → 200 Ω

### Ordering Information

Description	Type	Order number
SECULIFE   VL	Variable Load	M695V

For more information on accessories, see

- In the catalog *Measurement and Testing Technology*
- On the Internet at [www.gossenmetrawatt.com](http://www.gossenmetrawatt.com)